



Projector Stand

Hardware Documentation

1 April, 2026

Table of Contents

- 1 Product Overview 3
- 2 Product Details..... 4
- 3 Assembly and Installation 5
 - 3.1 What's included 5
 - 3.2 General recommendations 5
 - 3.3 Assembling the stand 5
- 4 Product Usage..... 6
- 5 Product Specifications 7
 - 5.1 General specifications 7
- 6 Cleaning and Maintenance..... 8
- 7 Troubleshooting your Device 9
- 8 Related Links 10
- 9 Compliance, Safety and Warranty Information 11

1 Product Overview

The VPixx Projector Stand is an adjustable-height stand designed to support the PROPixx projector or SHIELDPixx shielding enclosure in MRI and MEG environments. The stand is constructed entirely from MRI-safe, non-ferromagnetic materials. Features include:

- Adjustable height range: 30" to 43"
 - Constructed entirely from MRI-safe, non-ferromagnetic materials
 - Compatible with PROPixx projector and SHIELDPixx enclosure
-

2 Product Details

The stand consists of two vertical legs with height-adjustment collars and a top shelf that supports the PROPixx projector or SHIELDPixx enclosure.

Technical specifications:

- Adjustable height range: 30" to 43"
 - All materials are non-ferromagnetic and MRI-safe
-

3 Assembly and Installation

3.1 What's included

Qty	Item	MR safe?
1x	Projector stand (adjustable, MSR-safe materials)	Yes
1x	Allen key	Yes

3.2 General recommendations

Ensure that you observe the following before beginning assembly:

- Assemble this product outside the MSR. Never unpack and/or assemble your products in the room with the magnet.
- Two people are recommended when placing the SHIELDPixx on the assembled stand.

3.3 Assembling the stand

To adjust the top shelf height, use the Allen key to loosen the two screws on each leg collar. Slide the shelf to the desired height, and tighten the screws to lock the stand in place.

4 Product Usage

The projector stand supports the PROPixx projector or SHIELDPixx enclosure inside the magnetically shielded room. Once the stand is assembled and the projector or enclosure is in place, adjust the height to align the projection beam with the mirror and screen.

For complete installation instructions including projector placement, mirror alignment and image focusing, see the following product manuals:

- PROPixx Projector
 - SHIELDPixx Shielding Enclosure
 - Projection Mirror
-

5 Product Specifications

5.1 General specifications

- Stand height range: 30" to 43" (adjustable)
 - MR safety rating: MR Safe (non-ferromagnetic materials throughout)
 - Compatible with: PROPixx projector, SHIELDPixx shielding enclosure
-

6 Cleaning and Maintenance

The projector stand has a metal frame constructed from non-ferromagnetic materials. It can be cleaned with regular cleaning products or disinfectants. Spray the product into a soft cloth and wipe down the stand — do not directly spray the product on the stand.

7 Troubleshooting your Device

Issue	Solution
Stand is wobbly or unstable	Ensure all screws are fully tightened, including the height-adjustment collar screws.
Height adjustment is stiff	Loosen both screws on each leg collar before attempting to slide. Do not force the adjustment with only one screw loosened.

8 Related Links

- [PROPixx Projector](#)
 - [SHIELDPixx Shielding Enclosure](#)
 - [Projection Mirror](#)
 - [Projection Screens](#)
-

9 Compliance, Safety and Warranty Information

Click on the sections below to expand the relevant information.

Safety and Warnings

The following table shows MR safety warning indicators and their names and meanings. Please pay special attention to the symbols on your equipment, and obey all guidelines to protect staff and equipment from injury.

Symbol	Name	Description
	MR safe	Components marked with the green MR safe indicator present no danger to staff and equipment operating in the MR room.
	Conditionally MR safe	The yellow Conditionally MR safe indicator denotes components that are only MR safe in a limited fashion. Product-specific safety requirements apply. For additional details, consult the product manual's safety section.
	MR unsafe	The red MR unsafe indicator denotes components prohibited from entering the MR room.

Warranty Information

This product is warranted against manufacturing defects in materials and workmanship for two (2) years for parts and labour from the date of shipment.